

ISEO

This beautiful, delicate leather shows a luxurious look and elegant feel that reminds us of the aesthetic of calfskin. The smooth, slightly polished surface gains more subtle charm with age.

The Iseo embodies the essence of full aniline leather with higher visual appeal brought to a specific treatment that gives a translucent, polished surface. It features a palette inspired by the rich and warm earth tone colours that imparts a soft sense of elegance.

Custom Colors Available

CHARACTERISTICS

Size	4,5-5 mq / 50-55 sf
Thickness	1/1,3 mm
Raw material origin	European Bovine
Type of tannage	Chrome
Dyeing	Pure Aniline
Finishing	Aniline

PHYSICAL TEST

Rubbing resistance	
50 cycles (wet)	≥ 3 Grey Scale
20 cycles (sweat)	≥ 3 Grey Scale
Peeling strength	
	≥ 20 N/cm
Flex cracking	
	≥ 200.000 cycles (no cracking)
Light fastness	
	≥ 3 Blue Scale

LEATHER COMPLIANT TO

Cigarette test (BS 5852-1/79) – Match test
Cigarette test NFPA 260: 2003
IMO 652
IMVSS302
FS191 Method 5903 – Port Authority of NY & NJ
California TB 117 Section E

COMPLIANT ON REQUEST

FAR 25.853 (a)
FAR 25.853 (b)
UNI VF 9175 class1
IMO 653
BS 6853
DIN 5510
FS 191 Method 5903
Boston fire code BFD-IX-1

This leather is produced in conformity with REACH regulation

APPLICATIONS



AVAILABLE TREATMENTS

H2O 1
ACOUSTIC
FIRST CLASS



IS 01




IS 02



IS 03



IS 04




IS 05



SK 06



IS 07



IS 08



IS 09




IS 10



IS 12




IS 13



IS 14



IS 15



SK 16



SK 17



IS 18



IS 19



IS 20



IS 21



IS 22




IS 23



IS 24




IS 25




IS 26


AVAILABLE
EMBOSSED




A




B




C




D




E




F




G




H




I



L




M



N




O



P



Q



R



S



T



U



W



X

The standard placement of the plate are shown in the diagram below.


HIDE

Average size: 4,5-5mq /50-55 sf


HIDE

Plate size: 1000x1300 mm / 40x50 inch


OPTION 1




OPTION 2



OPTION 3



OPTION 4



The matchlines could be visible. Please verify the plate positioning prior to place an order.

MAINTENANCE ANILINE AND NUBUCK

- Avoid placing leather in direct sun
- Keep the leather away from direct heat
- Dust regularly with a soft and dry cloth
- Spots can be cleaned with a lightly damp cloth and mild soap
- Dry immediately after the treatment
- Do not rub. Use a light circular motion working toward the centre of the stain.

CLEANING ANILINE AND NUBUCK

The following guidelines will help you to maintain the beauty of nubuck or aniline without destroying any of the natural characteristics:

Dry carefully the area , and use one of the following methods.

Pre-testing on a hidden area.

Light cleaner

4Nubuck Foglizzo Light Cleaner will clean most stains and soiling from Nabuck leather. It is water based product that penetrates into the surface to clean without harming your Nabuck Leather.


Spray cream

4Nubuck Foglizzo Spray Cream will provide excellent protection for your Nubuck and Suede leathers against oil and water based spills.


Power Cleaner

4Nubuck Foglizzo Power Cleaner efficiently removes surplus fat and oil spots on Nubuck leather. It is to be used for oil based stains. Try first to clean the stain with 4Nubuck Foglizzo Cleaning Cloth and 4Nubuck Foglizzo Light Cleaner. If this is not successful, use this product.

PRODUCT FEATURES AND IMPERFECTIONS THAT MAY OCCUR



Neck creases



Backbone feature

Foglizzo Leather is a natural product, therefore every imperfections shown below or other hallmarks depends on its animal origins.
Please consider the images shown prior to placing an order.